

EMERGENCY GRAB BAG & DEFIBRILLATOR

Checklist weekly & after use

	Date				
Respirator					
Check the Powered Air Purifying Respirator (PAPR) battery is plugged in and hood and belt in the box. Ensure this is cleaned and decontaminated following use.					
Defibrillator & B.V.M					
Defibrillator Mindray D1 Check status indicator is Flashing Green. If red please contact EBME immediately					
Perform User Test (After Use)					
Check expiry date of battery (Monthly)					
1 set of defibrillator pads (in date) in sealed packet					
Bag – Valve – Mask					
RED pouch					
ECG Leads					
BLUE pouch					
ECG dots					
Tough Cut Scissors					
Disposable razor					
Grey zipped pocket					
1 set of defibrillator pads (in date) in sealed packet					
Tamper tags					
Top pocket					
Adult pocket mask					
15L reservoir mask					
Airway - Side pocket					
Oropharyngeal airway size 2, 3 & 4					
Nasopharyngeal airway size 6 & 7					
Sachet of lubricating gel					
Fluids - Side pocket					
1 x 500ml 0.9% Sodium Chloride					
1 x Wide Bore Giving Sets					
Back pocket					
Orange pouch					
Cannulae 4 (green, pink, blue)					
5 x Cannula dressings					
Green pouch					
6 x 5ml 0.9% Sodium Chloride (NaCl) (PosiFlush)					
Black pouch					
Tourniquet					
Chlorhexidine 2% wipes					
Gauze swabs					
1 x roll of surgical tape					
Front pocket					
3 x pairs gloves (medium)					
3 x pairs gloves (large)					
3 x pairs gloves (small)					
IGEL Size 3, 4 and 5					
Sharps bin					
Tamper Tag ID					
Nurse signature					