

**Paediatric Resus Trolley Checklist - including Daily top of trolley check & Weekly Checks and after every use**

Ward/Department **Childrens Unit** Month ... ..

**Daily Top of Trolley Check (incl. equipment on drip stand)**

Date	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Day (M,T,W, Th, F,S, Su)																															
Disposable gloves. Small medium and large																															
Small sharps bin																															
Resus trolley folder Including instructions for ligature cutter																															
Pocket mask 1x paediatric 1x adult																															
Stethoscope																															
Ambubags – 1 x child, infant and adult size with masks																															
O2 cylinder -green tubing attached with paediatric face mask																															
Anaesthetic breathing set																															
Surface clean and free from clutter																															
Portable suction (HDU only)																															
Initials																															



Paediatric Resus Trolley Checklist - including Daily top of trolley check & Weekly Checks and after every use

Daily Defibrillator Checklist

Date	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Day (M,T,W, Th, F,S, Su)																															
Mindray D3 Defibrillator plugged in, clean and dust free. No visible damage.																															
Battery installed with at least 60% of battery capacity																															
2 sets of defibrillator pads unopened in sealed packet (in date)																															
ECG leads present, undamaged, and connected securely.																															
ECG electrodes in date (1 packet)																															
Recorder paper present																															
Check Auto check completed and illuminated green tick visible.																															
If red cross, please complete User Test and report to EMBE if unresolved																															
Complete User Test: - Routine test - Energy deliver test																															
<b>Weekly &amp; after use</b>																															
Initials																															

Please file this form on your ward after completion and retain for minimum of 2 years

### 1<sup>st</sup> Drawer

<b>Number of removed security tag</b>					
<b>Date</b> Expiry dates to be checked on ALL items					
Laryngoscopes Miller 00 Please check the light works Mac 0 on each blade by attaching the Mac 1 handle. Mac 2 The packet does not need to Mac 3 be opened					
Oropharyngeal airways Sizes 000 x 1 00 x 1 0 x 1 1 x 2 1.5 x 2 2 x 2                    3 x 1                    4 x 1					
Nasopharyngeal airway size 4.0,4.5 5.0,5.5, 6.0 & 7.0					
Tongue depressor					
Endo tracheal (ET) tubes uncuffed sizes 2mm                    5mm 2.5mm                5.5mm 3mm                    6mm 3.5mm                6.5mm 4mm                    7 mm 4.5mm					
<b>Croup tube</b>					
CO <sub>2</sub> Monitor (HDU only)					
ET tube cuffed size 6mm, 7mm & 8mm					
Tracheal tube introducer s size 5, 4 and 2					
Angled catheter mount					
Aquagel x 5					
20ml syringe					
Tie for ET tube					
Elastoplast					
Duoderm					
Scissors					
Magill's forceps					
Yankuer suction catheters paediatric and adult 2 of each					
Ligature cutter					
Signature:					
Designation					

**Please file this form on your ward after completion and retain for minimum of 2 years**

**2<sup>nd</sup> Drawer**

Expiry dates to be checked on ALL items	Date	Date	Date	Date	Date
Paediatric and adult oxygen face mask					
Paediatric and adult nebuliser mask					
Paediatric oxygen non-rebreathing mask Adult oxygen rebreathing mask					
Enteral syringes 50ml , 20ml , 5 ml					
Ryles tube					
2 x Suction catheters size 5, 6, 8, 10,12,14					
Naso-gastric tubes sizes 6 x 2 8 x 2 10 x 2					
Oxygen saturation probe					
Signature					
Designation					

**Please file this form on your ward after completion and retain for minimum of 2 years**

### 3<sup>rd</sup> Drawer

Expiry dates to be checked on ALL items	Date	Date	Date	Date	Date
Tourniquet					
Blood bottles 2 x U and E, FBC, glucose coag., ESR 2x blood culture bottle 2 x Blood gas (capillary tubes)					
Razor					
Sterile stitch cutter and sterile blade					
4 grey ( 16g) cannulas      1x pink (20g) 4 green ( 18g) cannulas      1x white ( 14g) 4 safety cannulas size 25g 4 x yellow ( 24g) jelco 4 x blue (22g) jelco					
Safety Butterfly needles 2 x sizes 19g 23g and 25g					
2 x tenderfoot					
Hydrofilm/tegaderm x 4					
5 x 5ml Posiflushes					
Sterile gauze					
2 bandages					
Baby boards					
2 x 3 way taps 3 x extension sets 3 x double lumen extension sets					
Syringes 50ml x3 20ml x 4 10ml x 4 5ml x 4					

Expiry dates to be checked on ALL items	Date	Date	Date	Date	Date
2ml x 4 1ml x 4					
4 x filter needles 4 x blunt needles					
Capillary blood tubes					
10 x sanicloth					
Vaseline					
3 x chloraprep					
Drug additive labels					
4 x fluid free flow sets					
2 syringe pump giving sets					
Intraosseous gun, needle and dressings Check light is green					
Signature					
Designation					

**Please file this form on your ward after completion and retain for minimum of 2 years**

Paediatric Resus Trolley Combined Daily Checklist V4.1 2023 Reviewed October 2023 (inclusion of Defib checklist)

### 4th Drawer

Expiry dates to be checked on ALL items	Date	Date	Date	Date	Date
First line emergency drugs box BNF drug book Clock					
2 x fluid giving sets					
Anaphylaxis box					
1 x 500ml NaCl 0.45% with 5% dextrose 1 x 500ml NaCl 0.9% 1 x 500ml jelofusin 1 x 500ml dextrose 10% 1 x 500ml dextrose 5%					
NaCl 2.7% (HDU only)					
Portable suction (Ward only) Check unit is charged					
Security tag number used to lock the trolley					
Signature:					
Designation:					

**Please file this form on your ward after completion and retain for minimum of 2 years**